

AMPLIRUN® TOTAL RESPIRATORY VIRAL PANEL CONTROL (SWAB)

For *in vitro* diagnostic use

MBTC020: A panel of ten purified respiratory viruses pooled that have been inactivated to render them non-infectious and formulated in viral transport medium. Table 1 lists type of virus and cell-line used in the culture of each respiratory virus included in this control. This kit is intended to validate and control sample processing, amplification and detection in nucleic acid assays based on the molecular identification of respiratory viruses, using the product as an external run control.

INTRODUCTION:

Respiratory tract infections are a major cause of hospitalization. In humans these infections can be caused, in a great proportion, by a heterogeneous group of viruses that produce similar clinical presentations. Recent outbreaks associated to respiratory viruses highlight the importance of a rapid and accurate laboratory diagnosis.

CHARACTERISTICS:

The content is lyophilized. It is necessary to reconstitute it before use (refer to "Preparation of the reagents"). Total Controls are designed for single use, excess material should be discarded. Nucleic acid detection requires an extraction step that releases DNA/RNA for amplification and detection.

Product description:

Viral particles were purified from supernatants of infected cells by differential centrifugation (see Table 1 for cell-line used). Viruses were inactivated, rendering them non-infectious and diluted in viral transport medium containing cells obtained from epithelial human cell lines.

VIRUS	CELL LINE
ADENOVIRUS 4	HEp-2
CORONAVIRUS	MRC-5
INFLUENZA A H3N2	MDCK
INFLUENZA B	MDCK
NOVEL INFLUENZA A H1N1	MDCK
PARAINFLUENZA 1	LLC-MK2
PARAINFLUENZA 2	LLC-MK2
PARAINFLUENZA 3	LLC-MK2
RESPIRATORY SYNCYTIAL VIRUS (subtype A)	HEp-2
RESPIRATORY SYNCYTIAL VIRUS (subtype B)	HEp-2

Table 1.

KIT CONTENTS:

1 VIRCELL TOTAL RESPIRATORY VIRAL PANEL CONTROL (SWAB): 10 vials with a pooled of lyophilized respiratory viruses simulating a respiratory clinical sample. Each virus is in a concentration that ranges from 4000-10000 copies/vial. Batch concentration is provided in Certificate of Analysis.

Quantification validation was performed by real-time PCR.

Materials required but not supplied:

Molecular Biology grade water
Additional extraction and detection kit.

STORAGE REQUIREMENTS:

Special transport conditions not required. Store the lyophilized vial at 2-8°C. After reconstitution, suspension should be used on the same day. Unused product should be discarded.

STABILITY AND HANDLING OF REAGENTS:

Handle reagents in aseptic conditions to avoid microbial contaminations.

Use only the amount of reagent required for the test.

The product is stable until the expiry date indicated in the label, if the instructions for use are followed.

VIRCELL, S.L. does not accept responsibility for the mishandling of the reagents included in the kit.

RECOMMENDATIONS AND PRECAUTIONS:

1. This product is for *in vitro* diagnosis use only and for professional qualified staff.
2. Sterile tips with aerosol barrier are essential to prevent contamination.
3. Specimens should be handled as in the case of infectious samples using safety laboratory procedures. Thoroughly clean and disinfect all work surfaces with a freshly prepared solution of 0.5% sodium hypochlorite in deionized or distilled water.
4. In order to perform the test it is essential to have separate working areas.
5. Dispose of unused reagents and waste in accordance with all applicable regulations.
6. The component VIRCELL TOTAL CONTROL could include genetic material or substances of animal and/or human origin. VIRCELL TOTAL CONTROL contains inactivated microorganism, nevertheless, it should be considered potentially infectious and handled with care. Inactivation was verified by the absence of growth under same culture conditions used for each microorganism. No present method can offer complete assurance that these or other infectious agents are absent. All materials should be handled and disposed as of potentially infectious. Observe the local regulations for waste disposal.

PREPARATION OF THE REAGENTS:

1. Add 200 µl of Molecular Biology grade water to vial 1 and mix until completely reconstituted. The concentration will be approximately 35000 copies/ml for each virus in the panel once reconstituted.
2. Shake with vortex for 30 seconds to dissolve and homogenize completely.
3. Follow diagnostic kit instructions treating TOTAL CONTROL in an identical manner to a clinical specimen using recommended amount for extraction and detection.

INTERNAL QUALITY CONTROL:

Each batch is subjected to internal quality control testing before releasing. Quality control analysis is performed using a sample preparation kit and real-time PCR for quantification. Final quality control results for each particular lot are available.

INTERPRETATION OF RESULTS AND VALIDATION PROTOCOL FOR USERS:

Refer to indications of additional extraction and detection kit.

LIMITATIONS OF METHOD:

1. This reagent is intended to be used with methods of human diagnostics. This test has not been verified with other methods.
2. The user of this kit is advised to read carefully and understand the package insert. Strict adherence to the protocol is necessary to obtain reliable test results.
3. Use of this product should be limited only to personnel trained in molecular techniques.
4. This external run control does not substitute internal diagnostic kit controls.
5. Quantification conclusions cannot be drawn from a single point sample of known concentration. Precise clinical sample quantification could only be achieved by the standard curve method using a calibrator.
6. AMPLIRUN® TOTAL has not been designed to be used with a particular diagnostic kit coming from a certain manufacturer. It is used to validate and control sample processing, analysis and detection of a diagnostic laboratory functioning procedure.



PERFORMANCES:**• IDENTITY TEST**

PCR analysis was performed after extraction with a specific oligonucleotide pair for each virus present in the panel. The reactions produced fragments of the expected sizes.

• QUANTIFICATION TEST

Quantification is based on Real-Time qPCR using the standard curve method. This method involves the use of multiple replicates of different serial dilutions of both the product and the standard of quantification.










• INTRA-ASSAY PRECISION

3 vials of the product were extracted under identical extraction conditions and 3 replicas of each extraction were amplified by the same operator under identical qPCR conditions. Less than 15% coefficient of variance was observed between all assays.

• INTER-ASSAY PRECISION

1 vial of the product was extracted and 3 replicates from this vial were amplified by 2 different operators on 3 consecutive days. Less than 15% coefficient of variance was observed between all assays.

SYMBOLS USED IN LABELS:

	In vitro diagnostic medical device
	Use by (expiration date)
	Store at x-°C
	Batch code
	Catalogue number
	Consult instructions for use
	Reconstitute in x µl
	Shipment temperature
	Storage temperature

BIBLIOGRAPHY:

1. Bustin SA. (2000). Absolute quantification of mRNA using real-time reverse transcription polymerase chain reaction assays. *J Mol Endocrinol.* 25(2):169-193.
2. Bustin SA. (2002). Quantification of mRNA using real-time reverse transcription PCR (RT-PCR): trends and problems. *J Mol Endocrinol.* 29(1):23-39.
3. Chen Y1, Cui D, Zheng S, Yang S, Tong J, Yang D, Fan J, Zhang J, Lou B, Li X, Zhuge X, Ye B, Chen B, Mao W, Tan Y, Xu G, Chen Z, Chen N, Li L. (2009). Simultaneous detection of influenza A, influenza B, and respiratory syncytial viruses and subtyping of influenza A H3N2 virus and H1N1 virus by multiplex real-time PCR. *J Clin Microbiol.* 49(4):1653-6.
4. de-Paris F, Beck C, Machado AB, Paiva RM, da Silva Menezes D, de Souza Nunes L, Kuchenbecker R, Barth AL (2012). Optimization of one-step duplex real-time RT-PCR for detection of influenza and respiratory syncytial virus in nasopharyngeal aspirates. *J Virol Methods.* 186(1-2):189-92
5. Freeman WM, Walker SJ, Vrana KE. (1999). Quantitative RT-PCR: pitfalls and potential. *Biotechniques.* 1999 26(1):112-122, 124-125.
6. Gaunt ER, Hardie A, Claas EC, Simmonds P, Templeton KE (2010). Epidemiology and clinical presentations of the four human coronaviruses 229E, HKU1, NL63, and OC43 detected over 3 years using a novel multiplex real-time PCR method. *J Clin Microbiol.* Aug;48(8):2940-7.
7. Gunson RN1, Collins TC, Carman WF (2005). Real-time RT-PCR detection of 12 respiratory viral infections in four triplex reactions. *J Clin Virol.* 33(4):341-4.
8. Hartmann NM1, Dartscht M, Szewzyk R, Selinka HC (2013). Monitoring of adenovirus serotypes in environmental samples by combined PCR and melting point analyses. *Virol J.* Jun 10;10:190.
9. Larionov A, Krause A, Miller W. (2005). A standard curve based method for relative real time PCR data processing. *BMC Bioinformatics.* 6:62.
10. Osiowy C (1998). Direct detection of respiratory syncytial virus, parainfluenza virus, and adenovirus in clinical respiratory specimens by a multiplex reverse transcription-PCR assay. *J Clin Microbiol.* 36(11):3149-54.
11. Selvaraju SB1, Selvarangan R (2010). Evaluation of three influenza A and B real-time reverse transcription-PCR assays and a new 2009 H1N1 assay for detection of influenza viruses. *J Clin Microbiol.* 48(11):3870-5.
12. Syrmis MW, Whitley DM, Thomas M, Mackay IM, Williamson J, Siebert DJ, Nissen MD, Sloots TP (2004). A sensitive, specific, and cost-effective multiplex reverse transcriptase-PCR assay for the detection of seven common respiratory viruses in respiratory samples. *J Mol Diagn.* 6(2):125-31

For any question please contact:
customerservice@vircell.com

REVISED: 2019-01-15
L-MBTC020-EN-01

