AMPLIRUN® TOTAL RESPIRATORY VIRAL PANEL CONTROL (SWAB)

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BIBLIOGRAFIA:

1. Bustin SA. (2000). Absolute quantification of mRNA using real-time reverse transcription polymerase chain reaction assays. *J Mol Endocrinol.* 25(2):169-193.
2. Bustin SA. (2002). Quantification of mRNA using real-time reverse transcription PCR (RT-PCR): trends and problems. *J Mol Endocrinol.* 29(1):23-39.
3. Chen Y1, Cui D, Zheng S, Yang S, Tong J, Yang D, Fan J, Zhang J, Lou B, Li X, Zhuge X, Ye B, Chen B, Mao W, Tan Y, Xu G, Chen Z, Chen N, Li L. (2009). Simultaneous detection of influenza A, influenza B, and respiratory syncytial viruses and subtyping of influenza A H3N2 virus and H1N1 virus by multiplex real-time PCR. *J Clin Microbiol.* 49(4):1653-6.
4. de-Paris F, Beck C, Machado AB, Paiva RM, da Silva Menezes D, de Souza Nunes L, Kuchenbecker R, Barth AL (2012). Optimization of one-step duplex real-time RT-PCR for detection of influenza and respiratory syncytial virus in nasopharyngeal aspirates. *J Virol Methods.* 186(1-2):189-92
5. Freeman WM, Walker SJ, Vrana KE. (1999). Quantitative RT-PCR: pitfalls and potential. *Biotechniques.* 1999 26(1):112-122, 124-125.
6. Gaunt ER, Hardie A, Claas EC, Simmonds P, Templeton KE (2010). Epidemiology and clinical presentations of the four human coronaviruses 229E, HKU1, NL63, and OC43 detected over 3 years using a novel multiplex real-time PCR method. *J Clin Microbiol.* Aug;48(8):2940-7.
7. Gunson RN1, Collins TC, Carman WF (2005). Real-time RT-PCR detection of 12 respiratory viral infections in four triplex reactions. *J Clin Virol.* 33(4):341-4.
8. Hartmann NM1, Dartscht M, Szewzyk R, Selinka HC (2013). Monitoring of adenovirus serotypes in environmental samples by combined PCR and melting point analyses. *Virol J.* Jun 10;10:190.
9. Larionov A, Krause A, Miller W. (2005). A standard curve based method for relative real time PCR data processing. *BMC Bioinformatics.* 6:62.
10. Osioy C (1998). Direct detection of respiratory syncytial virus, parainfluenza virus, and adenovirus in clinical respiratory specimens by a multiplex reverse transcription-PCR assay. *J Clin Microbiol.* 36(11):3149-54.
11. Selvaraju SB1, Selvarangan R (2010). Evaluation of three influenza A and B real-time reverse transcription-PCR assays and a new 2009 H1N1 assay for detection of influenza viruses. *J Clin Microbiol.* 48(11):3870-5.
12. Syrmis MW, Whitley DM, Thomas M, Mackay IM, Williamson J, Siebert DJ, Nissen MD, Sloots TP (2004). A sensitive, specific, and cost-effective multiplex reverse transcriptase-PCR assay for the detection of seven common respiratory viruses in respiratory samples. *J Mol Diagn.* 6(2):125-31

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